



## Sterilizing and Disinfecting Products for Critical Environments

Sporicidin sterilizing and disinfecting products have been used in hospitals, medical offices, laboratories, and pharmaceutical and industrial cleanrooms since 1978. FDA cleared for disinfection of medical devices, Sporicidin disinfectants are also EPA registered for the disinfection of hospital environments and equipment, and compliant with the OSHA bloodborne pathogen standard.\*

For all of your sterilizing and disinfecting needs, turn to Sporicidin—a trusted name in hospitals and controlled environments for more than 30 years.

\* 29CFR 1910.1030



### Sporicidin® Disinfectant Solution and Towelettes

Cleans, disinfects and deodorizes

Sporicidin Disinfectant is an EPA registered intermediate level disinfectant that delivers a broad-spectrum kill that cleans, disinfects, and deodorizes. Sporicidin offers continuous bacteriostatic residual activity for up to 6 months. Alcohol-free, Sporicidin Brand Disinfectant is compatible with plastics, treated wood, glass, and metals. Sporicidin is non-staining, non-abrasive, and non-corrosive. Sporicidin provides users with a product safely rated in the lowest EPA toxicity category. It offers users the most neutral pH of any phenolic-based disinfectant. Sporicidin Disinfectant is EPA registered, OSHA compliant and available in ready-to-use gallon and spray bottles, as well as presaturated wipes.



Part No.	Description	Packaging
PS-3212F	Sporicidin Disinfectant, 32 oz. (0.9L), with trigger spray, fresh scent	12 bottles/case
RE-1284F	Sporicidin Disinfectant, 1 gallon (3.8L) container, fresh scent	4 gallons/case
CAN-18012C	Sporicidin Wipes, 5" x 8" (12.7 x 20.3cm), citrus scent	180 wipes/canister; 12 canisters/case
JUM-8506C	Sporicidin Wipes, 9.5" x 12" (24.1 x 30.5cm), citrus scent	85 wipes/canister; 6 canisters/case
ST50-600	Sporicidin Wipes, 5" x 8" (12.7 x 20.3cm), individually wrapped, citrus scent	50 wipes/box; 12 boxes/case

### Sporicidin® Antimicrobial Soap

Alcohol-free lotion soap

The Sporicidin Antimicrobial Soap has been categorized by the FDA as a Surgical Scrub and Health Care Personnel hand wash. It is ideal for ISO class 6 and 7 buffer zones, as well as ISO Class 8 Anterooms found in hospitals and compounding pharmacies. The soap kills up to 99.9% of gram positive and gram-negative pathogenic and odor-causing organisms. Alcohol-free Sporicidin Soap can be used daily for hand washing prior to donning gloves or beginning work in controlled environments or isolators.



Part No.	Description	Packaging
CALS-2416	Sporicidin Antimicrobial Soap, 16 oz. (0.5L), with pump	24 bottles/case
ALS-1282	Sporicidin Antimicrobial Soap, 1 gallon (3.8L) container	2 gallons/case



P.O. Box 530 • Spartanburg, SC 29304 USA  
 Tel: +1-864-503-8333  
 Toll Free (US, Canada, PR): 1-800-424-3733  
 Europe: +33 297437690  
[www.sporicidin.com](http://www.sporicidin.com) • [sporicidin@contecinc.com](mailto:sporicidin@contecinc.com)



**A trusted name in Laboratories and Controlled Environments**

*Sporicidin disinfectants have been used in hospitals, medical offices, laboratories, and pharmaceutical and industrial clean rooms since 1978.*

**Sporicidin® Sterilizing and High Level Disinfecting Solution**

Sporicidin Sterilizing and High Level Disinfecting Solution is intended to sterilize or provide high level disinfection for reusable medical instruments and equipment. An alkaline glutaraldehyde activator combined with a phenol/phenate buffer solution, it is effective with the lowest level of glutaraldehyde cleared by the FDA to support efficacy claims (1%). With a 3-year shelf life (non activated) and 14-day active use life, this disinfecting solution is a cost effective option for medical device sterilization.



Part No.	Description	Packaging
SSS-1284C	Sporicidin Sterilizing Solution, 1 gallon (3.8L) container	4 gallons/case
SSS-3216C	Sporicidin Sterilizing Solution, 32 oz. (0.95L) container	4 bottles/box; 4 boxes/case
PMG-6000	Sporicidin Indicator Kit; 60 glutaraldehyde strips & 60 phenol strips	4 kits/case

**Sporicidin® Enzymatic Cleaner**

Sporicidin Enzymatic Cleaner is a concentrated solution with a blend of enzymes that effectively cleans instruments and surfaces soiled with organic and inorganic debris. This highly effective, low-foaming cleaner and stain remover has a neutral pH, is biodegradable, non-flammable, and non-toxic. It is ideal for use in ultrasonic cleaners and reduces time spent on gross filth removal. One gallon makes 128 gallons of ready-to-use cleaner.



Part No.	Description	Packaging
ENZ-1284	Sporicidin Enzymatic Cleaner, 1 gallon (3.8L) container	4 gallons/case

**Sporicidin® Gamma Irradiated Disinfectant Solution**

Sporicidin Gamma Irradiated Disinfectant is double-bagged and gamma irradiated for pharmaceutical, biotech, and medical device clean room operations. The product is sterile lot released according to the USP <71> protocol and provides continuous bacteriostatic residual activity for up to 6 months.



Part No.	Description	Packaging
GAMMA-1284C	Sporicidin Gamma Irradiated Disinfectant Solution, 1 gallon (3.8L) container	4 gallons/case

**Ask about our new Contec® Hand Hygiene System:**

- Contec® Foam Hand Sanitizer
- Contec® Alcohol-free Foam Hand Sanitizer



P.O. Box 530 • Spartanburg, SC 29304 USA  
 Tel: +1-864-503-8333  
 Toll Free (US, Canada, PR): 1-800-424-3733  
 Europe: +33 297437690  
[www.sporicidin.com](http://www.sporicidin.com) • [sporicidin@contecinc.com](mailto:sporicidin@contecinc.com)