

## Sporicidin Material Compatibility

<b>Plastics</b>	Polypropelene Polyurethane Teflon Polyacetal Polysulfonate Polyetherimide Polycarbonate Vinyl chloride Polybutylene terephthalate (PET – Polyester) Polypphenylene oxide
<b>Synthetic rubbers</b>	Silicone rubber Acrylonitrile-butadiene rubber (NBR) Neoprene Nitrile rubber Fluoro rubber
<b>Glass</b>	Optical glass
<b>Alumite</b>	Anodized Aluminum
<b>Stainless Steel</b>	Yes
<b>Adhesive</b>	Epoxy adhesive
<b>Sealing Material</b>	Silicone sealing compound
<b>Dental Impression Materials</b>	Royalite ABS and ABS/PVC plastic extruded sheets Irreversible Hydrocolloid Silicone Condensation Silicone Vinyl Polysiloxane Jeltrate Alginate Full Arch Alginate Type IV Gypsum Polysulfide Rubber Base Polyether